

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-235
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. WAVEX – CLASSIC TUBE ROTATOR

Valid for Catalog No. E11300

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description:

- 1.1.1 Wavex tube rotators provide gentle but effective mixing of biological samples, held in 1.5 ml to 50 ml tubes.
- 1.1.2 Used in a variety of applications including immune precipitations, prevention of blood coagulation & latex diagnostics.
- 1.1.3 Wavex-Classic Tube Rotator has Disk Type Design.

1.2 Design:

Disk type design

1.3 performance Specifications:

1.3.1 Voltage (AC)	: 220~240 Volt
1.3.2 Frequency (Hz)	: 50
1.3.3 Power (W)	: 20
1.3.4 Tilt Angle	: 0- 90°
1.3.5 Motor type	: DC Motor
1.3.6 Speed range (RPM)	: 0 ~ 80
1.3.7 Operation time	: Continuous
1.3.8 Dimensions	: 300x220x310mm
1.3.9 Permissible ambient Temp	: 5- 40°

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-235
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. WAVEX – CLASSIC TUBE ROTATOR

1.4.0 Permissible Relative humidity: 80 %

1.4.1 Speed Display : Scale

1.4 Material: Body – ABS, Stainless Steel

1.5 Color: Dark Blue

1.6. Spare Accessories:

E11301 Plastic clips for 1.5ml x 12 MCT, Disk accessory – 12pcs/case
--

E11302 Plastic clips for 15ml x 12 MCT, Disk accessory – 12pcs/case

E11303 Plastic clips for 50ml x 12 MCT, Disk accessory – 12pcs/case

1.7 CE certified

1.8 Chemical resistance: see website (www.abdoslifesciences.com)

1.9 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing

Instruction for Use

Storage temp. Max 30°C

Before using the product check crackness, etc.

Check all the information printed on pouches or boxes.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-235
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
WAVEX – CLASSIC TUBE ROTATOR

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

Manufactured & Marketed By :
Abdos Labtech Private Limited

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com
Lic. No.: MFG/IVD/2024/000115 MADE IN INDIA



Scan to visit website

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110



TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-235
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
WAVEX – CLASSIC TUBE ROTATOR

