

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-245
Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
WAVEMAX – Digital Incubator

Valid for Catalog No. E11355, E11356

SRN (Single Registration Number): IN-MF-000022272

1. Description Specifications:

1.1 Description: WAVEMAX – Digital Incubator

1.1.1 Wavemax-Digital Incubators provide the facility to monitor the chamber temperatures in real time.

1.1.2 Has large internal chambers with two adjustable shelves for increased space. Internal Power Outlet for Installation of other

1.1.3 Equipment for sample agitation within the Digital Incubators.

1.1.5 Wavemax-Digital Incubator is available in a heating unit or heating and cooling for lower temperature applications.

1.2 Design: Digital Incubator

1.3 performance Specifications:

Specifications	E11355	E11356
Voltage	220-240V	220-240V
Frequency	50-60 Hz	50-60 Hz
Temp. range	Ambient +5 to 60°C	Ambient +5 to 60°C
Platform Dimensions	235x290mm	235x290mm
Operating Temp. range	+4 to 65°C	+4 to 65°C
Exterior	335x370x475mm	335x370x475mm
Interior	26x23.5x32.5 cm	26x23.5x32.5 cm
Capacity	20L	20L

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-245
Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. **WAVEMAX – Digital Incubator**

Weight	6.5kg	6.5kg
Temp. Increment	1°C	1°C
Temp. Operation	Only Heating	Heating & Cooling

1.4 Material: Body – ABS

1.5 Color: Grey

1.6 E11355 – Wavemax – Heating Digital Incubator

1.7 E11355 – Wavemax – Heating & Cooling Digital Incubator

1.8 Spare Accessories

1.8.1 **E11357** – Extra shelf 10.5x8 inch plastic (not autoclavable)

1.8.2. CE certified

1.8.3 Chemical **resistance:** see website (www.abdoslifesciences.com)

1.9 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing.

Instruction for Use

Storage temp. Max 30°C

Before using the product check crackness, etc.

Check all the information printed on pouches or boxes.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-245
Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
WAVEMAX – Digital Incubator

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

Manufactured & Marketed By :
Abdos Labtech Private Limited

Plot No. 1 Lakshewari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com
Lic. No.: MFG/IVD/2024/000115 MADE IN INDIA

IVD

UK
CA



Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

