

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-163  
Rev No: 04 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Dual Purpose Scoop, PS**

**Valid for Catalog No. P11907, P11908**

**SRN (Single Registration Number): IN-MF-000022272**

## 1. Description/Specification

### 1.1 Description: Dual Purpose Scoop, PS

1.1.1 ABDOS Dual purpose Scoop is designed to combine the function of a scoop, weigh boat, and funnel to make this product economically beneficial.

1.1.2 These Scoops are made from high quality anti-static polystyrene that is waterproof and reduces the probability of materials sticking to the scoop.

1.1.3 These scoops have locking tabs that help to retain materials.

1.1.4 These scoops can withstand powders, gels, and semi-solids.

1.1.5 The broad tip design keeps the scoop stable while on the scale and enables materials to be easily scooped and dispensed.

### 1.2 Design: Round bottom Scoop

### 1.3 Performance Specifications:

Code	Size	Dimension (mm)	Colour	No. Per Case
P11907	Large	127x90	White	100
P11908	Small	127x63	White	100

### 1.4 Material:

1.4.1 made up of Hydrophobic Polystyrene (PS), Confirm with US FDA 21CFR

1.4.2 Free from Heavy metals, it is safe to use in contact with food stuff, Pharmaceuticals & drinking water.

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-163  
Rev No: 04 Date- 05/2024-25

## **Abdos Labtech Pvt. Ltd. Dual Purpose Scoop, PS**

---

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

**1.5 Color:** Translucent (Natural)

**1.6 Quality control:** Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release amendments.

### **1.7 Other information**

1.7.1 **Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

1.7.2 **Manufactured** in Aseptic Environment

1.7.3 **Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

**1.8 Sterilization:** NA

### **Instruction for Use**

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

### **WARNING**

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.


**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-163  
Rev No: 04 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Dual Purpose Scoop, PS**



Quality Manager  
Abdos Labtech Private Limited  
Employee Code 30110

