

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-168
Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Funnel Holder, PP

Valid for Catalog No. P70311, P70312

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Funnel Holder, PP

1.1.1 These are corrosion free & can hold funnels with diameter 3 inc to 6 inches.

1.1.2 Extra plate provided can be placed to hold funnels with diameter 1 to 3 inches.

1.1.3 Two tapered wedges are provided to have firm grip on rods.

1.2 Design: Funnel Holder, PP

1.3 Performance Specifications:

Code	Places	No. per Case
P70311	Single	6
P70312	Double	6

1.4 Material:

1.4.1 PP- Polypropylene (PP), confirm US FDA 21CFR

1.4.2 Free from Heavy metals, it is safe to use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-168
Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Funnel Holder, PP

1.5 Color: White

1.6 Quality control: Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release amendments.

1.7 Basic Features

1.7.1 Made up of High pure virgin Polypropylene (PP)

1.7.2 Temperature range: -80⁰C to +121⁰C

1.7.3 Autoclavability: Yes, Autoclavable at 15 psi and 121⁰C for 20 minutes

1.7.4 Chemical resistance: See website website (www.abdoslifesciences.com)

1.7.5 Manufactured in Aseptic Environment

1.7.6 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Autoclaving Guideline for Plastics

All plastic consumables are autoclavable at 121⁰C for 20 minutes, 15 psi/1atm. Do not use a dry cycle as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to at least 50⁰C before removing product.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30⁰C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-168
Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Funnel Holder, PP

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

UKRP
Obelis UK LTD
Sandford Gate,
East Point Business Park,
Oxford OX4 6LB
United Kingdom
Tel: +44 (0)1865 378012
info@obelis.co.uk

CE **IVD** **UKCA** **REACH** **RoHS**

Obelis s.a.
Bd General Wafsis 53
1030 Brussels Belgium
Tel: +32 (0) 2 732 55 54
Fax: +32 (0) 2 732 60 03
mail@obelis.net

NON STERILE **30°C**

Manufactured & Marketed By:
Abdos Labtech Private Limited
Plot No. 1, Lakheshwari Bhagwanpur Roorkee-247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91 011-41081897 or Email us at labtech@abdosindia.com
Lic. No.: MFG/IVD/2024/000155 MADE IN INDIA

Electronic II U
Scan to visit website

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110



TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-168
Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Funnel Holder, PP

