

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-173
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Assorted Workstation Storage Bin, Acrylic

Valid for Catalog No. P70104

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Assorted Workstation Storage Bin, Acrylic.

1.1.1 Assorted Workstation Storage Bin can organize most of their commonly used laboratory disposables in one location.

1.1.2 First section of bin is used to store tubes, vials and pipette tips.

1.1.3 Second section of the bin is used to store and dispense tissue paper, Kim wipes and Para film M (supplied with super-sharp serrated blade)

1.1.4 Third section of the bins is used to store and dispense gloves.

1.2 Design: Workstation Storage Bin

1.3 Performance Specifications:

Code	No. of Compartment	L x W x H (mm)	No. Per Case
P70104	5	255 x 135 x 368	1

1.4 Material:

1.4.1 Poly (methyl methacrylate) (PMMA) Acrylic

1.4.2 Free from Heavy metals, It is safe use in contact with food stuff, Pharmaceuticals & drinking water

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-173
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Assorted Workstation Storage Bin, Acrylic

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color: Natural

1.6 Quality control: Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release amendments.

1.7 Other information:

1.7.1 Strong and lightweight material.

1.7.2 Good impact strength, higher than both glass & polystyrene.

1.7.3 Transmits up to 92% of visible light (3 mm thickness)

1.7.4 Better chemical stability.

1.8. Basic Features

1.8.1 Autoclavability: NO, Non autoclavable

1.8.2 Chemical resistance: see website (www.abdoslifesciences.com)

1.8.3 Manufactured in Aseptic Environment

1.7.6 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-173
Rev No: 07 Date- 05/2024-25


Abdos Labtech Pvt. Ltd. **Assorted Workstation Storage Bin, Acrylic**

Instruction for Use

Product can be used for food contact also.
Storage temp. Max 30°C
Before using the product check crackness etc.
Check all the information printed on pouches or boxes.
Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure
Consult local fire codes prior to storage of flammable liquids.



Quality Manager

Abdos Labtech Private Limited

Employee Code 30110



TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-173
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Assorted Workstation Storage Bin, Acrylic

