

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-175  
Rev No: 07 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Drying Rack, SS**

**Valid for Catalog No. P20302**

**SRN (Single Registration Number):** IN-MF-000022272

## 1. Description/Specification

### 1.1 Description: Drying Rack

1.1.1 These Drying Racks are used for draining & drying test tubes, beakers, flasks and other lab ware.

1.1.2 There are 55 pegs in the rack, 50 smaller pegs for products with small diameter & 5 large bracket designed to hold Erlenmeyer flasks, beakers etc.

1.1.3 The pegs are plastic coated for better protection.

1.1.4 The rack can be placed on the rack, on the table or near the sink.

1.1.5 The rack comes with a tray at the bottom on order to avoid drops on the surface or table.

### 1.2 Design: Drying Rack with 55 pegs

### 1.3 Performance Specifications:

Code	Places	L x W x H (mm)	No. Per Case
P20302	55	360 x 160 x 680	1

### 1.4 Material: Epoxy Coated Steel

### 1.5 Color: White

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-175  
Rev No: 07 Date- 05/2024-25

## **Abdos Labtech Pvt. Ltd. Drying Rack, SS**

---

**1.6 Quality control:** Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release amendments.

### **1.7 Other information**

1.7.1 Strong and lightweight material.

1.7.2 Good impact strength

1.7.3 Better chemical resistance.

### **1.8. Basic Features**

**1.8.1 Autoclavability:** autoclavable, but these racks are used for drying purpose.

**1.8.2 Chemical resistance:** see website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.8.3 Manufactured** in Aseptic Environment

**1.8.4 Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

### **Instruction for Use**

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

### **WARNING**

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.


**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-175  
Rev No: 07 Date- 05/2024-25

## Abdos Labtech Pvt. Ltd. Drying Rack, SS



Quality Manager  
Abdos Labtech Private Limited  
Employee Code 30110

