

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-194

Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Magnetic Retrievers, PTFE/LNICO V Magnet

Valid for Catalog No.

P80187, P80188, P80189

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Magnetic Retrievers, PTFE/LNICO V Magnet

1.2 Design: Magnetic Retrievers

1.3 Performance Specification:

1.3.1 Used to retrieve stirrer bar.

1.3.2 Inert & leak proof

1.3.3 Strong alloy magnets sheathed in pure white inert PTFE.

1.4 Material:

1.4.1 Moulded in from FDA approved PTFE & LNICO v Magnet.

1.5 Color: White

1.6 Quality Control: Testing of various characteristics in accordance with the valid Specification.

1.7 Basic Features:

1.7.1 Made up of Highly Pure Virgin PTFE Material.

1.7.2 Resistance to virtual all chemical attacks.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-194

Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Magnetic Retrievers, PTFE/LNICO V Magnet

Code	Type	ODx L (mm)	No./ Case	Basic Features
P80187	Magnetic Retrievers, PTFE/LNICO V Magnet	10x150	2	- Used to retrieve stirrer bar.
P80188		10x250	2	-Resistance with chemicals & with stand upto 260 ⁰ C
P80189		10x350	2	

1.8 Basic Features:

1.8.1 Chemical resistance: see website (www.abdoslifesciences.com)

1.8.2 Manufactured in Aseptic Environment

1.8.3 Autoclavability: Yes

1.8.4 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-194

Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Magnetic Retrievers, PTFE/LNICO V Magnet

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

