

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-36
Rev No: 01 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Starter Packs-Elegant & Premium Plus

Valid for Catalog No. E11148, E11149, E11159 E11160

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Starter Elegant & Premium Plus Pipettes, PP

1.1.1 Pipettes are fully autoclavable.

1.1.2 Pipettes are lightweight, ergonomic, with low force design.

1.1.3 Volume settings can be clearly read through digital display.

1.1.4 Pipettes are easy to calibrate & maintain with the tool supplied.

1.1.5 Pipettes are designed to avoid repetitive strain injuries.

1.1.6 Pipettes are supplied with operator manual that provides complete information on calibration & maintenance of the pipettes.

1.1.7 Pipettes are suitable for transfer liquid volumes.

1.1.8 Matching universal style pipette tips, calibrated in accordance with ISO 8655. The quality control involves gravimetric testing of each pipette with distilled water at 22°C.

1.2 Design: Light weight & simple to use operation.

1.3 Performance Specifications: Push Button, Cap, Push Button Rod, Tip ejector, Grip cover, Handle, Display cover, Tip ejector collar, Tip cone, Volume setting, Sealing & ejecting tips.

1.4 Material: PP (polypropylene)

1.5 Color: White

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-36
Rev No: 01 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Starter Packs-Elegant & Premium Plus

1.6 Basic Features

1.6.1 Made up of highly Virgin thermoplastic polypropylene.

1.6.2 Abdos premium plus pipettes are fully autoclavable.

1.6.3 CE certified

1.6.4 Autoclavability: Yes

1.6.5 Chemical resistance: see website (www.abdoslifesciences.com)

1.7. Pipettes

Code	Description	Capacity	No. Per Case
E11148	Premium Plus Starter Pack 01	Premium Plus Starter Pack 01, Includes: E11107 (0.1-10µl Pipette), E11109 (10-100µl Pipette), E11111 (200-1000µl Pipette), P10139 (0.1-10µl Gen2 Tip Box), P10140 (1-200µl Gen2 Tip Box), P10109 (100-1000µl Gen2 Tip Box), P70201 (5 Place Linear Pipette Stand), Pipette Bench-top Hangers, Abdos Biodegradable Pen & Pencil Set	1
E11149	Premium Plus Starter Pack 02	Premium Plus Starter Pack 02, Includes: E11108 (2-20µl Pipette), E11110 (20-200µl Pipette), E11111 (200-1000µl Pipette), P10139 (0.1-10µl Gen2 Tip Box), P10140 (1-200µl Gen2 Tip Box), P10109 (100-1000µl Gen2 Tip Box), P70201 (5 Place Linear Pipette Stand), Pipette Bench-top Hangers, Abdos Biodegradable Pen & Pencil Set	1
E11159	Elegant Starter Pack 01	Elegant Starter Pack 01, Includes: E11151 (0.1-10µl Pipette), E11153 (10-100µl Pipette), E11156 (100-1000µl Pipette), P10139 (0.1-10µl Gen2 Tip Box), P10140 (1-200µl Gen2 Tip Box), P10109 (100-1000µl Gen2 Tip Box), P70201 (5 Place Linear Pipette Stand), Pipette Benchtop Hangers, Abdos Biodegradable Pen & Pencil Set	1
E11160	Elegant Starter Pack 02	Elegant Starter Pack 02, Includes: E11152 (2-20µl Pipette), E11154 (20-200µl Pipette), E11156 (100-1000µl Pipette), P10139 (0.1-10µl Gen2 Tip Box), P10140 (1-200µl Gen2 Tip Box),	1

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-36
Rev No: 01 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Starter Packs-Elegant & Premium Plus

		P10109 (100-1000µl Gen2Tip Box), P70201 (5 Place Linear Pipette Stand), Pipette Benchtop Hangers, Abdos Biodegradable Pen & Pencil Set	
--	--	--	--

1.8 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing

Instruction for Use

Product can be used for food contact also.
Storage temp. Max 30°C
Before using the product check crackness etc.
Check all the information printed on pouches or boxes.
Always wear gloves while handling specimen.

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure
Consult local fire codes prior to storage of flammable liquids.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-36
Rev No: 01 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Starter Packs-Elegant & Premium Plus

Manufactured & Marketed By :
Abdos Labtech Private Limited

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com
Lic. No.: MFG/IVD/2024/000115 MADE IN INDIA



Scan to Visit website

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110



TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-36
Rev No: 01 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Starter Packs-Elegant & Premium Plus

