

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-40
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Pipette Rack Stand, Acrylic

Valid for Catalog No. P70201

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Pipettes Rack Stand, Acrylic

1.1.1 Keep your pipettes within easy inside this convenient rack.

1.1.2 Holds five standard size pipettes.

1.1.3 Racks are mounted on their own stand.

1.1.4 Keep working area or desk free from clutter.

1.2 Design: Pipette Rack Stand

1.3 Performance Specifications:

Code	Places	L x W x H (mm)	No. per Case
P70201	5	190x135x265	1

1.4 Material:

1.4.1 Acrylic, Poly (methyl methacrylate) PMMA, confirm US FDA 21CFR

1.4.2 Free from Heavy metals, It is safe use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.3 Product is manufactured without the use of Raw Materials of animal origin Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color: Transparent

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-40
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Pipette Rack Stand, Acrylic

1.6 Quality control: Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7. Basic Features

1.7.1 Made up of Acrylic, Poly (methyl methacrylate) PMMA

1.7.2 Strong & Lightweight material, Plate Thickness approx. 7.5mm

1.7.3 Good impact strength, higher than both glass and polystyrene.

1.7.4 Fitment with ABDOS Pipettes only

1.7.5 **Autoclavability:** NO

1.7.6 Manufactured in aseptic environment.

1.7.7 **Chemical resistance:** see website (www.abdoslifesciences.com)

1.8 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-40
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Pipette Rack Stand, Acrylic

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure.

Consult local fire codes prior to storage of flammable liquids.

Manufactured & Marketed By:

Abdos Labtech Private Limited

Plot No. 1, Lakhewari, Bhagwanpur, Roorkee-247667 (UK) India

For Complaints/Feedback Contact Customer Care at above address or:

Ph: +91-011-41081897 or Email us at: labtech@abdosindia.com

Lic. No.: MFG/IVD/2024/000158 MADE IN INDIA



Scan to visit website

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

