

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-43
Rev No: 09 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Pipette Tips, Bulk PP

Valid for Cat. No. P10115, P10137, P10129, P10130, P10101, P10105, P10138, P10102, P10106, P10174, P10175

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Disposable PP Micropipette Tips, DNase, RNase, Human DNA, ATP & Pyrogen Free, CE & IVD certified.

1.2 Design:

1.2.1 Open Cylindrical Upper and conical Apex with orifice for liquid dispensing

1.3 Performance Specifications:

Code	Capacity	Colour	No. per Case
P10115	0.1-10µl (Graduated)	Natural	1000
P10137	0.1-10µl Extended (Graduated)	Natural	1000
P10129	1-200µl (Bevelled Graduated)	Yellow	1000
P10130	1-200µl (Bevelled Graduated)	Natural	1000
P10101	1-200µl (Xtreme Point Tip Orifice)	Yellow	1000
P10105	1-200µl (Xtreme Point Tip Orifice)	Natural	1000
P10138	1-300µl (Graduated)	Natural	1000
P10102	100 - 1000µl	Blue	500
P10106	100 - 1000µl	Natural	500
P10174	100 - 1000µl XL	Natural	500

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-43
Rev No: 09 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Pipette Tips, Bulk PP

P10175	100 - 1250µl	Natural	500
---------------	--------------	---------	-----

1.3.2 Fitting to pipettes: Universal Fit (For detail refer the Abdos compatibility chart)

1.4 Material:

1.4.1 Polypropylene (PP), Hydrophobic, confirm US FDA 21CFR

1.4.2 The product complies with USP class VI. Free from Heavy metals, it is safe to use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color:

Tip color- Natural, Yellow, and Blue

1.6 Sterilization: NA

1.7 Quality control: Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.8 Other information: Recommended for single use only

1.9 Basic Features:

1.9.1. Made up of: High Purity Virgin Polypropylene (PP)

1.9.2. Temperature range: -80⁰C to + 121⁰C for 20 minutes

1.9.3. Autoclavability: Yes, Autoclavable at 15 psi and 121 ⁰C for 20 minutes

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-43
Rev No: 09 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Pipette Tips, Bulk PP

Quality Manager
Abdos Labtech Private Limited
Employee Code 30110



P10115



P10137



P10129



P10130



P10138



P10102



P10175

