

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 71
Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Handypette Pipette Aid.

Valid for Catalog No. E11218, E11216, E11222

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Handypette Pipette Aid.

1.1.1 Handypette Pipette Aid ensures quick & accurate operation.

1.1.2 Colour Coded by volume with green for 10 ml, red for 25 ml & blue for 2ml.

1.1.3 Users can easily insert glass & plastic pipette into the collar with stabilizing fingers for secure hold.

1.1.4 Resistant to acids, alkalis & easily disassembly for cleaning.

1.1.5 Rotating the knurled thumbwheel draws liquid to desired volumes. To release liquid, apply light pressure to the fast release lever. For graduation dispensing rotate the thumb wheel.

1.2 Design: Handypette Pipette

1.3 Performance Specifications:

Code	Capacity (ml)	No. per Case
E11218	10	4
E11216	25	4
E11222	2	4

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 71
Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Handypette Pipette Aid.

1.4 Quality control: Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.5. Basic Features:

1.5.1 Autoclavability: NO, Non autoclavable

1.5.2 Chemical resistance: see website (www.abdoslifesciences.com)

1.5.3 CE certified

1.5.4 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 71
Rev No: 06 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Handypette Pipette Aid.



Quality Manager
Abdos Labtech Private Limited
Employee Code 30110

