

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 72  
Rev No: 07 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Rubber Pipette Pump**

**Valid for Catalog No. E11217**

**SRN (Single Registration Number):** IN-MF-000022272

## 1. Description/Specification

### 1.1 Description: Rubber Pipette Pump

- 1.1.1 Rubber Pipette Pump ensures quick & accurate operation.
- 1.1.2 Pump opening is standard to fit almost all glass & plastic pipette.
- 1.1.3 One size can be used for liquid up to 100 ml.
- 1.1.4 The pump has a heavy wall and easy to grip body.

### 1.2 Design: Rubber Pipette Pump

### 1.3 Performance Specifications:

| Code   | Capacity (ml) | No. per Case |
|--------|---------------|--------------|
| E11217 | Up to 100ml   | 4            |

### 1.4 Quality control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

### 1.5. Basic Features:

**1.5.1 Autoclavability:** NO, Non autoclavable

**1.5.2 Chemical resistance:** see website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 72  
Rev No: 07 Date- 05/2024-25

## Abdos Labtech Pvt. Ltd. Rubber Pipette Pump

**1.5.3 Colour-** Available in various colour


**1.5.4 Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

### Instruction for Use

Product can be used for food contact also.  
Storage temp. Max 30°C  
Before using the product check crackness etc.  
Check all the information printed on pouches or boxes.  
Always wear gloves while handling specimens.

### WARNING

**DO NOT** store strong oxidizing agents.  
**DO NOT** place any plastic labware in a flame.  
**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure  
Consult local fire codes prior to storage of flammable liquids.



Quality Manager  
Abdos Labtech Private Limited  
Employee Code 30110



# TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 72  
Rev No: 07 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Rubber Pipette Pump**

---

