

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 73
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Rave Pipette Controller, PP

Valid for Catalog No. E11113

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Rave Pipette controller ideal for Repetitive Pipetting performance
Maximizing Pipetting comfort

1.2 Design: Robust, Lightweight & simple to use operation.

1.3 Performance Specifications:

1.3.1 Can use volumetric or serological pipettes ranging from 0.1 ml to 100 ml.

1.3.2 It contains integral & replacement 0.3 µm Filter.

1.3.3 Autoclavable silicon pipette holder

1.4 Material : PP (polypropylene), confirm US FDA 21CFR

1.5 Color : Grey

1.6 Basic Fetures

1.6.1 Made up of highly Virgin thermoplastic polypropylene

1.6.2 Autoclavability: Yes (Only Silicon Pipette Holder – Don't autoclave whole instrument)

1.6.3 CE certified

1.6.4 Chemical resistance: See website (www.abdoslifesciences.com)

1.6.5 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 73
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Rave Pipette Controller, PP

1.7 Dimension:

Code	Dimension
E11113	: Length- 210 mm, ϕ -60mm
E11119	: (Spare) 3.0 μ m Filter

1.8 Quality control: Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing

Instruction for Use

Product can be used for food contact also.
Storage temp. Max 30°C
Before using the product check crackness, pouch punctured etc.
Check all the information printed on pouches or boxes.
Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure
Consult local fire codes prior to storage of flammable liquids.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 73
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Rave Pipette Controller, PP

Manufactured & Marketed By :
Abdos Labtech Private Limited

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com
Lic. No.: MFG/IVD/2024/000115 MADE IN INDIA

IVD

UK
CA

CE

Electronic IFU



Scan to visit website

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

