

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 77  
Rev No: 06 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Pipette Storage Rack with Magnet, ABS**

**Valid for Catalog No. P70202**

**SRN (Single Registration Number): IN-MF-000022272**

## 1. Description

1.1.1. Pipette Storage Rack with Magnet, ABS

1.1.2 These racks help to protect & organize glass and plastic pipettes.

1.1.3 Angled four shelved compartments keep 0.1 to 50 ml pipettes or cans within easy reach while preventing damage and clutter.

1.1.4 These have cut-out viewing window and come with two removable attachments that hold pipette bulbs, fillers and pumps.

1.1.5 These unique innovative magnetized racks have strong magnets on both sides and mounts securely to most metal surfaces.

**1.2 Design:** Pipette Storage Rack with Magnet, ABS

**1.3 Performance Specifications:**

Code	Description	L x W x H (mm)	No. per Case
<b>P70202</b>	Storage 1-50ml pipette or cans	291x88x401	1

## 1.4 Material:

1.4.1 Acrylonitrile Butadiene Styrene (ABS), confirm US FDA 21CFR

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 77  
Rev No: 06 Date- 05/2024-25

## **Abdos Labtech Pvt. Ltd.** **Pipette Storage Rack with Magnet, ABS**

---

1.4.2 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

### **1.5 Color:** Blue

**1.6 Quality control:** Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

### **1.6. Basic Features**

**1.6.1 Autoclavability:** NO, Non autoclavable

**1.6.2 Chemical resistance:** see website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.6.3 Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

### **Instruction for Use**

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 77  
Rev No: 06 Date- 05/2024-25

## Abdos Labtech Pvt. Ltd. Pipette Storage Rack with Magnet, ABS


### WARNING

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.

**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure

Wear a face shield and / or safety goggles if you have to work anywhere near a centrifuge that's in Consult local fire codes prior to storage of flammable liquids.



Quality Manager  
Abdos Labtech Private Limited  
Employee Code 30110

