

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 78
Rev No: 04 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Pipette Storage Rack, ABS

Valid for Catalog No. P10620

SRN (Single Registration Number): IN-MF-000022272

1. Description

1.1.1. Pipette Storage Rack, ABS

1.1.2 These racks can be assembled horizontally or vertically to build the desired configuration, or used as individual holders, using sections.

1.1.3 Angle shelve design makes these racks compatible to use for storing 0.1 to 50ml of serological pipettes.

1.1.4 Can be used bench top or can be connected to most metal surfaces.

1.2 Design: Pipette Storage Rack, ABS

1.3 Performance Specifications:

Code	Compatible	W x H x D (mm)	No. per Case
P10620	0.1to 50ml	301x100x127.5	3

1.4 Material:

1.4.1 Acrylonitrile Butadiene Styrene (ABS), confirm US FDA 21CFR

1.4.2 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color: Blue

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 78
Rev No: 04 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Pipette Storage Rack, ABS

1.6 Quality control: Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.6. Basic Features

1.6.1 Autoclavability: NO, Non autoclavable

1.6.2 Chemical resistance: see website (www.abdoslifesciences.com)

1.6.3 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Product can be used for food contact also.
Storage temp. Max 30°C
Before using the product check crackness etc.
Check all the information printed on pouches or boxes.
Always wear gloves while handling specimens.







WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure
Consult local fire codes prior to storage of flammable liquids.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH- 78
Rev No: 04 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Pipette Storage Rack, ABS

					UKRP Obelis UK LTD Surreyford Gate, East Point Business Park, Oxford OX4 6LB United Kingdom Tel: +44 1491 378012 info@obelis.co.uk
Manufactured & Marketed By: Abdos Labtech Private Limited Plot No. 1, Lakshewari Bhagwanpur Roorkee-247667 (UK) India For Complaints/Feedback Contact Customer Care at above address or Ph: +91-011 41081897 or Email us at labtech@abdosindia.com Lic. No.: MFG/IVD/2024/000158 MADE IN INDIA					 Scan to visit website

Quality Manager
Abdos Labtech Private Limited
Employee Code 30110

