

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ Micro-83- II

Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Petri Dish, PS- Gama Radiated Sterile

Valid for Cat. No. P10905, P10906, P10907, P10908, P10909, P10910, P10935, P10936, P10937, P10938

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Petri Dish- Gama Radiated Sterile, DNase, RNase, Human DNA, ATP, Pyrogen Free, CE & IVD certified, ISO24998

1.2 Design: Shallow dish design

1.3 Performance Specifications:

Code	Diameter (mm)	Height (mm)	Compartment	Vented	Packing	No. Per Case
P10905	90	15.80	No	3 Vent	Sleeve of 20	480
P10906	90	15.80	No	3 Vent	Sleeve of 10	480
P10907	90	15.80	No	3 Vent	Individual Wrapped	380
P10908	90	15.80	No	0 Vent	Sleeve of 20	480
P10909	90	15.80	No	0 Vent	Sleeve of 10	480
P10910	90	15.80	No	0 Vent	Individual Wrapped	380
P10935	60	14.20	No	0 Vent	Sleeve of 20	480
P10936	60	14.20	No	3 Vent	Sleeve of 20	480

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ Micro-83- II

Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Petri Dish, PS- Gama Radiated Sterile

P10937	140	17.50	No	0 Vent	Sleeve of 5	180
P10938	140	17.50	No	3 Vent	Sleeve of 5	180

1.4 Material:

1.4.1 Petri Dish made up of Hydrophobic Polystyrene (PS), Confirm with US FDA 21CFR

1.4.2 The product complies with USP class VI. Free from Heavy metals, it is safe to use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color: Transparent (Natural color)

1.6 Quality control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7 Other information

1.7.1 Petri Dishes are suitable for all bacteriological applications.

1.7.2 Vented dishes have elevated lid that allows for superior gas exchange and evaporation.

1.7.3 Petri dishes retain their shape up to 55°C.

1.7.4 Easy to Stack, 12° inclination 10 stacking filed with 20 ml media for 90mm.

1.7.5 Recommended for single use only

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ Micro-83- II

Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Petri Dish, PS- Gama Radiated Sterile

1.8 Basic Features

1.8.1 Temperature range: Stable up to 60°C.

1.8.2 Autoclavability: NO

1.8.3 Chemical resistance: See website (www.abdoslifesciences.com)

1.8.4 Manufactured in Clean room ISO Class 8 (100,000)

1.8.5 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

1.9 Sterilization: Sterile catalog number from this lot have been pre-sterilizes and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products- Requirements for validation and routine control- Radiation sterilization*. The product has been sterilized with Co-60 Gamma Radiation. SAL 10^{-6}

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness, pouch punctured etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ Micro-83- II

Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Petri Dish, PS- Gama Radiated Sterile

    		UKRP Obelis UK LTD Sandford Gate, East Point Business Park, Oxford OX4 6LB United Kingdom Tel: +44 1491 378012 info@obelis.co.uk
Manufactured & Marketed By : Abdos Labtech Private Limited Plot No. 1 Lakshwari Bhagwanpur Roorkee-247667 (UK) India For Complaints/Feedback Contact Customer Care at above address or Ph: +91-011-41081897 or Email us at: labtech@abdosindia.com Lic. No.: MFG/IVD/2024/000115 MADP IN INDIA		Obelis s.a. Bel General Walsis S3 2030 Brussels Belgium Tel: +32 2 732 55 54 Fax: +32 2 732 60 63 mail@obelis.net  Electronic I U Scan to visit website

Quality Manager 
Abdos Labtech Private Limited
Employee Code 30110


