

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-197
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Resealable Bags with Zip Lock, PE

Valid for Catalog No. U40110, U40111

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Resealable Bags with Zip Lock, PE

1.1.1 These bags can be used to store and ship samples, solids and food made of 52 micron/210 gauge thickness.

1.1.2 Suitable for use with food.

1.1.3 Though it can be used for all applications it has especially used in food, water, forensic science, and quality control laboratory.

1.1.4 Using the rim and channel closure the bag can be closed by figure pressure and reopened repeatedly.

1.1.5 Content stored in the bag would be dust free and free of contamination.

1.2 Design:

Rectangular shaped Resealable Bags.

1.3 Performance Specifications:

Code	Volume (Ltr)	Size WxL (inch)	Weight (gm)	No. Per Case
U40110	4	9.8 x 10.2	6	100
U40111	7	13 x 13.35	11	100

1.4 Material:

1.4.1 PE (Polyethylene), Hydrophobic, confirm US FDA 21CFR

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-197
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Resealable Bags with Zip Lock, PE

1.4.2 Free from Heavy metals, It is safe use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color: Transparent Natural

1.6 Quality control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7 Other information

1.7.1 Better chemical stability.

1.7.3 Highly transparent, elastic thermoplastic.

1.7.4 **Autoclave** - NO

1.7.5 **Chemical resistance:** See website (www.abdoslifesciences.com)

1.7.6 Manufactured in Aseptic Environment

1.7.7 **Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-197
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Resealable Bags with Zip Lock, PE

Instruction for Use

Product can be used for food contact also.
Storage temp. Max 30°C
Before using the product check crackness, pouch punctured etc.
Check all the information printed on pouches or boxes.
Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure.
Consult local fire codes prior to storage of flammable liquids.

Manufactured & Marketed By :

Abdos Labtech Private Limited

Plot No. 1 Lakshewari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com
Lic. No.: MFG/MD/2024/000898 MADE IN INDIA



Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

