

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-201  
Rev No: 07 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Safety Face Shield, PC**

**Valid for Catalog No. U40102**

**SRN (Single Registration Number): IN-MF-000022272**

## 1. Description/Specification

### 1.1 Description: Safety Face Shield, PC

1.1.1 Clear face shield of polycarbonate provides high impact and splash protection.

1.1.2 Impact resistant and lightweight frames for comfortable Face safety.

1.1.3 Ratched adjustment on the headgear provides comfortable fit.

### 1.2 Design: Safety Face Shield

### 1.3 Performance Specifications:

Code	L x W x H (cm)	No. Per Case
<b>U40102</b>	28.5 x 22.7 x 13	1

### 1.4 Material: PC (Polycarbonate)

**1.5 Color: Top: Black Frame: Transparent Natural**

### 1.6 Quality Control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

### 1.7 Other Information

1.7.1 Excellent clarity and visibility.

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-201  
Rev No: 07 Date- 05/2024-25

## **Abdos Labtech Pvt. Ltd. Safety Face Shield, PC**

1.7.2 Moderate chemical stability.

1.7.3 Transmits light better than alternative material.

**1.8 Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.9** Manufactured in Aseptic Environment

**2.0 Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

### **Instruction for Use**

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

### **WARNING**

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.

**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

**Marketed By :**

**Abdos Labtech Private Limited**

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India

For Complaints/Feedback Contact Customer Care at above address or

Ph.: +91-011-41081897 or Email us at [labtech@abdosindia.com](mailto:labtech@abdosindia.com)



Electronic IFU



Scan to visit website

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-201  
Rev No: 07 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Safety Face Shield, PC**

---

Quality Manager  
Abdos Labtech Private Limited  
Employee Code 30110

