

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-203  
Rev No: 09 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Abdos Derma Gaurd Blue Nitrile Gloves – 9.5 inches length**

**Valid for Catalog No. U20341, U20342, U20343**

**SRN (Single Registration Number): IN-MF-000022272**

## 1. Description/Specification

### 1.1 Description: Abdos Brand – Derma Gaurd Blue Nitrile Gloves – 9.5 inches length

1.1.1 These gloves are ideal products to protect your exposure to bacteria, viruses and chemicals. It also enhances your ability to offer maximum security against cross contamination.

1.1.2 Its unique blue color provides quick, visual differentiation from natural rubber latex and other synthetic gloves.

1.1.3 Ambidextrous, textured fingertips and beaded cuff enhances wet & dry grip.

1.1.4 Custom design enhances comfort & fit.

1.1.5 These gloves are powder free, minimizing the potential for powder related complications.

1.1.6 Conforms to ASTM D6319, D6124 and EN 455 Standards, EU standard directive 93/42/EEC ,Class I, IMDR 2027 (MD)

1.1.7 Derma Gaurd Blue Nitrile Gloves Powder free- sample meets the overall migration requirement under Commission Regulation (EU) No 10/2011- Plastic materials and articles shall not transfer their constituents to foodstuff in quantities exceeding 10mg of total constituents released per dm<sup>2</sup> of food contact surface (mg/dm<sup>2</sup>) (overall migration limit)

1.1.8 Palm Thickness :0.08mm (3.15mil) & Finger Thickness: 0.11mm (4.33mil)

1.1.9 Average Length: 240mm

### 1.2 Design: Hand gloves design

### 1.3 Performance Specifications: 4mil

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-203  
Rev No: 09 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Abdos Derma Gaurd Blue Nitrile Gloves – 9.5 inches length**

Code	Size	No. Per Case
U20341	Small	100
U20342	Medium	100
U20343	Large	100

## 6mil

Code	Size	No. Per Case
U20355	Small	100
U20356	Medium	100
U20357	Large	100

**1.4 Material:** Nitrile gloves

**1.5 Color:** Blue

**1.6 Quality control:**

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release

**1.7 Other information**

1.7.1 These gloves are used for working on lab, machine and any other work.

1.7.2 These gloves are ideal products to protect you from any of microbial contamination.

1.7.3 These Flexible gloves are fit in hand & easy to use on working.

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-203  
Rev No: 09 Date- 05/2024-25

## **Abdos Labtech Pvt. Ltd.** **Abdos Derma Gaurd Blue Nitrile Gloves – 9.5 inches length**

- 1.7.4 Powder free Chlorinated gloves
- 1.7.5 Suitable for food handling
- 1.7.6 **Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))
- 1.7.7 Manufactured in Aseptic Environment
- 1.7.8 **Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

### **Instruction for Use**

Product can be used for food contact also.  
Storage temp. Max 30°C  
Before using the product check crackness, punctured etc.  
Check all the information printed on pouches or boxes.  
Always wear gloves while handling specimens.

### **WARNING**

**DO NOT** store strong oxidizing agents.  
**DO NOT** place any plastic labware in a flame.  
**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure  
Consult local fire codes prior to storage of flammable liquids.

#### **Manufactured & Marketed By :**

#### **Abdos Labtech Private Limited**

Plot No. 1 Lakeshwari Bhagwanpur Roorkee 247667 (UK) India  
For Complaints/Feedback Contact Customer Care at above address or  
Ph.: +91-011-41081897 or Email us at [labtech@abdosindia.com](mailto:labtech@abdosindia.com)  
Lic. No.: MFG/MD/2024/000898 MADE IN INDIA



# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-203  
Rev No: 09 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Abdos Derma Gaurd Blue Nitrile Gloves – 9.5 inches length**

Quality Manager  
Abdos Labtech Private Limited  
Employee Code 30110

