

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-204
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Abdos Derma Gaurd White Latex Gloves – 9.5 inches length

Valid for Catalog No. U20344, U20345, U20346

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Abdos Brand – Derma Guard White Latex Gloves – 9.5 inches length

1.1.1 These gloves are ideal products to protect your exposure to bacteria, viruses and chemicals. It also enhances your ability to offer maximum security against cross contamination.

1.1.2 Ambidextrous, textured fingertips and beaded cuff enhances wet & dry grip.

1.1.3 Softness provides superior comfort & natural fit.

1.1.4 These gloves are powder free, minimizing the potential for powder related complications, irritation, and dermatitis.

1.1.5 Lower level of extractable proteins and chemical residue, hence lower incidences of allergic reactions.

1.1.6 Conforms to ASTM D3578 and EN 455 Standards, EU standard directive 93/42/EEC Class I

1.1.7 Palm Thickness :0.10mm (3.93mil) & Finger Thickness: 0.10mm (3.93mil)

1.1.8 Average Length: 240mm

1.2 Design: Hand gloves design.

1.3 Performance Specifications:

Code	Size	No. Per Case
U20344	Small	100

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-204
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Abdos Derma Gaurd White Latex Gloves – 9.5 inches length

U20345	Medium	100
U20346	Large	100

1.4 Material: Latex gloves

1.5 Color: White

1.6 Quality control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release

1.7 Other information

1.7.1 These gloves are used for working on lab, machine, and any other work.

1.7.2 These gloves are ideal products to protect you any of microbial contamination.

1.7.3 These Flexible gloves are fit in hand & easy to use on working.

1.7.4 Powder free Chlorinated gloves

1.7.5 **Chemical resistance:** See website (www.abdoslifesciences.com)

1.7.6 Manufactured in Aseptic Environment

1.7.7 **Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness, punctured etc.

Check all the information printed on pouches or boxes.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-204
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Abdos Derma Gaurd White Latex Gloves – 9.5 inches length

Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

Manufactured & Marketed By :

Abdos Labtech Private Limited

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com
Lic. No. : MFG/MD/2024/000898 MADE IN INDIA



Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

