

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-205
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
KIMTECH Purple/Blue Nitrile Gloves 9.5- & 12.0-inches length

Valid for Catalog No.

U20301, U20301B, U20302, U20302B, U20303, U20303B, U20304, U20305, U20306

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Kimtech Purple/Blue Nitrile Gloves

1.1.1 These gloves are ideal products to protect you exposure to bacteria, viruses and chemicals. It also enhances your ability to offer maximum security against cross contamination.

1.1.2 Ambidextrous, textured fingertips and beaded cuff enhances wet & dry grip.

1.1.3 Softness provides superior comfort & natural fit.

1.1.4 These gloves are powder free, minimizing the potential for powder related complications.

1.1.5 Kimtech Blue Nitrile Gloves Powder free- sample meets the overall migration requirement under Commission Regulation (EU) No 10/2011- Plastic materials and articles shall not transfer their constituents to foodstuff in quantities exceeding 10mg of total constituents released per dm² of food contact surface (mg/dm²) (overall migration limit)

1.1.6 Conforms to ASTM D6319 and EN 455 Standards, EU standard directive 93/42/EEC Class I

1.1.7 Palm Thickness :0.12mm (4.7mil) & Finger Thickness: 0.15mm (5.9mil)

1.1.8 Average Length: 240mm

1.2 Design: Hand gloves design

1.3 Performance Specifications:

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-205
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
KIMTECH Purple/Blue Nitrile Gloves 9.5- & 12.0-inches length

1.3.1 Kimtech Purple Nitrile Gloves – 9.5 inches length

Code	Size	No. Per Case
U20301	Small	100
U20302	Medium	100
U20303	Large	100

1.3.2 Kimtech Purple Nitrile Gloves – 12.0 inches length

Code	Size	No. Per Case
U20304	Small	50
U20305	Medium	50
U20306	Large	50

1.3.3 Kimtech Blue Nitrile Gloves – 9.5 inches length

Code	Size	No. Per Case
U20301B	Small	100
U20302B	Medium	100
U20303B	Large	100

1.4 Material: Nitrile gloves

1.5 Color: Purple & Blue

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-205
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. KIMTECH Purple/Blue Nitrile Gloves 9.5- & 12.0-inches length

1.6 Quality Control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7 Other information

1.7.1 These gloves are used for working on lab, machine and any other work.

1.7.2 These gloves are ideal products to protect you any of microbial contamination.

1.7.3 These Flexible gloves are fit in hand & easy to use on working.

1.7.4 Powder free Chlorinated gloves

1.7.5 Suitable for food handling

1.7.6 **Chemical resistance:** See website (www.abdoslifesciences.com)

1.7.7 Manufactured in Aseptic Environment

1.7.8 **Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness, punctured etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-205
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
KIMTECH Purple/Blue Nitrile Gloves 9.5- & 12.0-inches length

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

Marketed By :

Abdos Labtech Private Limited

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com



Scan to visit website

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

