

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-206
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Glove Box Holder, Acrylic

Valid for Catalog No. P70301, P70302, P70303

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Glove Box Holder, Acrylic

1.1.1 They designate a place for gloves.

1.1.2 Clear and pristine holder to give clean appearance.

1.1.3 Top loading of the glove box.

1.1.4 Easy to dispense the gloves through the holder.

1.1.5 Each holder has holes that help it to be installed in a convenient location to the proximity of your work and task.

1.2 Design: Glove Box Holder

1.3 Performance Specifications:

Code	No. of Places	L x W x H (mm)	No. Per Case
P70301	1	260 x 100 x 138	1
P70302	2	260 x 100 x 270	1
P70303	3	260 x 100 x 415	1

1.4 Material: Acrylic

1.5 Color: Transparent Natural

1.6 Quality control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-206
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Glove Box Holder, Acrylic

Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7 Other information

1.7.1 Excellent clarity and visibility.

1.7.2 Strong and lightweight material.

1.7.3 Good impact strength.

1.7.4 **Chemical resistance:** See website (www.abdoslifesciences.com)

1.7.5 Manufactured in Aseptic Environment

1.7.6 **Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.


DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-206
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Glove Box Holder, Acrylic



Quality Manager
Abdos Labtech Private Limited
Employee Code 30110

ABDOS LABTECH PVT. LTD.
QUALITY
DEPT.
ROORKEE, INDIA