

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-211
Rev No: 09 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Parafilm M Dispenser, Acrylic

Valid for Catalog No. P70401

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description:

Parafilm M dispenser – Dispensing 2 “and 4” Rolls.

Ideal for storing, dispensing, and cutting Parafilm® M tapes and labels.

1.2 Design: Rectangular Box size with easy open lid

1.3 Performance Specifications:

1.3.1 Dispenser is designed to store, dispense, and cut Parafilm®.

1.3.2 The Parafilm® dispenser has unique design for the cutting mechanism providing maximum safety.

1.4 Material:

1.4.1 Made up of Acrylic, confirm US FDA 21CFR.

1.4.2 Free from Heavy metals, It is safe use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color: Light grey

1.6 Quality control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-211
Rev No: 09 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Parafilm M Dispenser, Acrylic

Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7. Basic Features:

1.7.1 Made up of Highly Virgin Material.

1.7.2 The safety knife has unique design for cutting the film.

1.7.3 Dimensions: Length x width x Height (mm)
180 x 120 x 160mm

1.7.4 Capacity: 2 & 4 inches

1.7.5 Chemical resistance: See website (www.abdoslifesciences.com)

1.7.6 Manufactured in Aseptic Environment

1.7.7 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Storage temp. Max 30°C

Before using the product check crackness, pouch punctured etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-211
Rev No: 09 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Parafilm M Dispenser, Acrylic



Manufactured & Marketed By :
Abdos Labtech Private Limited

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com
Lic. No. : MFG/IVD/2024/000115 MADE IN INDIA



Scan to visit website

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

