

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-212  
Rev No: 07 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Parafilm M Dispenser, ABS**

**Valid for Catalog No. P70402**

**SRN (Single Registration Number): IN-MF-000022272**

## 1. Description/Specification

### 1.1 Description:

Parafilm M dispenser – Dispensing 2 “and 4” Rolls

Ideal for storing, dispensing, and cutting Parafilm® M tapes and labels.

**1.2 Design:** Rectangular Box size with easy open lid

### 1.3 Performance Specifications:

**1.3.1** Dispenser is designed to store, dispense, and cut Parafilm®.

**1.3.2** The Parafilm® dispenser has unique design for the cutting mechanism providing maximum safety.

### 1.4 Material:

**1.4.1** Acrylonitrile butadiene styrene (ABS) confirm US FDA 21CFR

**1.4.2** Free from Heavy metals, It is safe use in contact with food stuff, Pharmaceuticals & drinking water.

**1.4.3** Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

**1.5 Color:** Blue

### 1.6 Quality Control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-212  
Rev No: 07 Date- 05/2024-25

## **Abdos Labtech Pvt. Ltd. Parafilm M Dispenser, ABS**

Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

### **1.7. Basic Features:**

**1.7.1** Made up of Highly Virgin Material.

**1.7.2** The safety knife has unique design for cutting the film.

**1.7.3 Dimensions:** Length x Width x Height (mm)  
156 x 120 x 156 mm

**1.7.4 Capacity:** 2 & 4 inches

**1.7.5 Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.7.6** Manufactured in Aseptic Environment

**1.7.7 Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

### **Instruction for Use**

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

### **WARNING**

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.

**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure.

Consult local fire codes prior to storage of flammable liquids.

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-212  
Rev No: 07 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Parafilm M Dispenser, ABS**



**Marketed By :**  
**Abdos Labtech Private Limited**  
Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India  
For Complaints/Feedback Contact Customer Care at above address or  
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com



Quality Manager  
Abdos Labtech Private Limited  
Employee Code 30110

