

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-213
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Freezer Thermometer

Valid for Catalog No. U20610

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Freezer Thermometer

1.1.1 These thermometers can be used as probe to investigate and protect frozen and chilled samples from potential spoilage.

1.1.2 Monitor temperatures inside your freezer or refrigerators with the added convenience of sound & flashing LED lights to alert out of range temperatures.

1.1.3 Suction cup attaches the probe inside freezer/refrigerator while the thermometer secures to the outside with hook & loop tape.

1.1.4 Large digital display shows freezer/refrigerator temperature plus ambient temperature.

1.2 Design: Digital Thermometer

1.3 Performance Specifications:

Code	L x W x H (mm)	No. Per Case
U20610	136 x 19 x 28	1

1.4 Color: White

1.5 Quality Control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-213
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Freezer Thermometer

1.6 Other information:

1.6.1 Set freezer range: (-9.5° to -23°C).

1.6.2 Set refrigerator range: (2° to 8°C).

1.6.3 Can be used as Fahrenheit °F or Celsius °C

1.6.4 Includes one AAA battery.

1.6.5 Probe cable measures 0.9 m; probe measures 5 mm OD x 17 mm L

1.6.6 **Chemical resistance:** See website (www.abdoslifesciences.com)

1.6.7 Manufactured in Aseptic Environment

1.6.8 **Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-213
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Freezer Thermometer

Marketed By :
Abdos Labtech Private Limited
Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com



Quality Manager
Abdos Labtech Private Limited
Employee Code 30110

