

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-220  
Rev No: 07 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Multipurpose Labeling Tape**

**Valid for Catalog No.**

**U10701 W, U10701 B, U10701 G, U10701O, U10701 Y, U10701 R**

**SRN (Single Registration Number): IN-MF-000022272**

## 1. Description/Specification

### 1.1 Description: Multipurpose labeling Tape

1.1.1 These autoclavable tapes adhere to almost any dry surface-glass, plastic, metal, paper and rubber.

1.1.2 These tapes stick throughout the most demanding procedures, are oil-proof, and can be steam autoclaved.

1.1.3 Holds in a wide range of temperatures and transfers cleanly from one surface to another, leaving no messy residue.

1.1.4 These labels accept handwritten notations inscribed with pen, pencil or marker.

### 1.2 Design: Multipurpose labeling Tape Roll

### 1.3 Performance Specifications:

Code	Size (Inch)	Colour	No. Per Case
<b>U10701 W</b>	0.75 x 500	White	1 Roll
<b>U10701 B</b>	0.75 x 500	Blue	1 Roll
<b>U10701 G</b>	0.75 x 500	Green	1 Roll
<b>U10701 O</b>	0.75 x 500	Orange	1 Roll
<b>U10701 Y</b>	0.75 x 500	Yellow	1 Roll
<b>U10701 R</b>	0.75 x 500	Yellow	1 Roll

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-220  
Rev No: 07 Date- 05/2024-25

## **Abdos Labtech Pvt. Ltd. Multipurpose Labeling Tape**

---

**1.4 Color:** White, Blue, Green, Orange, Yellow, Red

### **1.5 Quality Control:**

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

### **1.6 Other information:**

**1.6.1 Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.6.2** Manufactured in Aseptic Environment

**1.6.3 Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

### **Instruction for Use**

Storage temp. Max 30°C  
Before using the product check crackness etc.  
Check all the information printed on pouches or boxes.  
Always wear gloves while handling specimens.

### **WARNING**

**DO NOT** store strong oxidizing agents.  
**DO NOT** place any plastic labware in a flame.  
**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure.  
Consult local fire codes prior to storage of flammable liquids.

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-220  
Rev No: 07 Date- 05/2024-25

## Abdos Labtech Pvt. Ltd. Multipurpose Labeling Tape

**Marketed By :**

**Abdos Labtech Private Limited**

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India  
For Complaints/Feedback Contact Customer Care at above address or  
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com



Electronic IFU



Scan to visit website

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

