

ABDOS

LIFE SCIENCES

Advancing Science with Quality



Universal Pipette Tips

SOLUTION TO YOUR PIPETTING NEEDS



Features

- ✧ Manufactured from ultrapure USP Class VI certified virgin polypropylene conforming to USFDA 21 CFR
- ✧ Range: 10-10,000 μ L
- ✧ Full traceability
- ✧ Colour coded racked tip box for visual identification
- ✧ Precise diamond polished bevelled tips prevents sample hold up at tip orifice
- ✧ Xtreme point tips for precise gel loading & careful pipetting between space as 1 mm in PAGE | AGE | IEF
- ✧ Packaging: Bulk | Individually Wrapped | Refill | Racked Tip Box
- ✧ Inner radial rings for perfect air lock grip with the pipettes
- ✧ Colour coded wafers in racked tip box have details that include product code, product & volume description with lot number. Colour coded wafers keep your work organised by visual identification

Bulk

- Easy-to-open and reseal
- Fully autoclavable
- Available in non-sterile



Individually Wrapped

- Fully automated blister packaging
- Gamma irradiated with SAL 10⁻⁶
- Complete traceability with mention of product code, lot number and expiry date



Racked Tip Box

- Available sterile and non-sterile
- Colour coded rack inserts
- Design patented generation 2 tip box for stability
- Full traceability




Refill

- Saves precious time and storage space
- Refill inserts can be easily changed



Quality & Compliance



FREE



LEAK TEST



FREE



FREE



FREE



FREE

SAL 10⁻⁶



GC/MS
Tested

IVDR
CERTIFIED

STERILE R

Ordering Information

Product Code	Product Description	No. Per Case
P10115	Micro Tips, Bulk in Resealable Bags, Capacity, 0.1-10 µl (Graduated) Natural	1000
P10137	Micro Tips, Bulk in Resealable Bags, Capacity, 0.1-10 µl Extended (Graduated) Natural	1000
P10129	Micro Tips, Bulk in Resealable Bags, Capacity, 1-200 µl (Bevelled Graduated) Yellow	1000
P10130	Micro Tips, Bulk in Resealable Bags, Capacity, 1-200 µl (Bevelled Graduated) Natural	1000
P10101	Micro Tips, Bulk in Resealable Bags, Capacity, 1-200 µl (Xtreme Point Tip Orifice) Yellow	1000
P10105	Micro Tips, Bulk in Resealable Bags, Capacity, 1-200 µl (Xtreme Point Tip Orifice) Natural	1000
P10138	Micro Tips, Bulk in Resealable Bags, Capacity, 1-300 µl (Graduated) Natural	1000
P10102	Micro Tips, Bulk in Resealable Bags, Capacity, 100-1000 µl (Graduated) Blue	500
P10106	Micro Tips, Bulk in Resealable Bags, Capacity, 100-1000 µl (Graduated) Natural	500
P10174	Micro Tips, Bulk in Resealable Bags, Capacity, 100-1000µl XL (Extended) (Graduated)	500
P10175	Micro Tips, Bulk in Resealable Bags, Capacity, 100-1250 µl (Graduated)	500
P10144	Macro Tips, Bulk in Resealable Bags, Capacity, ml 1-5 ml (Graduated) Natural	100
P10145	Macro Tips, Bulk in Resealable Bags, Capacity, ml 1-10 ml (Graduated) Natural	100
P10164	PureRight Individually Wrapped, Sterile, 0.2-10 µl (Graduated)	350
P10165	PureRight Individually Wrapped, Sterile, 0.2-10 µl Extended (Graduated)	350
P10161	PureRight Individually Wrapped, Sterile, 2-200 µl (Bevelled Graduated)	350
P10166	PureRight Individually Wrapped, Sterile, 2-300 µl (Graduated)	350
P10163	PureRight Individually Wrapped, Sterile, 100-1000 µl	350
P10167	PureRight Individually Wrapped, Sterile, 100-1000 µl XL (Graduated)	100
P10168	PureRight Individually Wrapped, Sterile, 100-1250 µl (Graduated)	100
P10169	PureRight Individually Wrapped, Sterile, 100-5000 µl (Graduated)	60
P10170	PureRight Individually Wrapped, Sterile, 500-10000 µl (Graduated)	50
P10116	Racked Micro Tips 0.2-10 µl (Graduated) Natural, No. Per Case -10 Racks x 96 Pcs	960
P10139	Racked Micro Tips 0.2-10 µl Extended (Graduated) Natural, No. Per Case -10 Racks x 96 Pcs	960
P10140	Racked Micro Tips 2-200 µl (Bevelled Graduated) Yellow, No. Per Case -10 Racks x 96 Pcs	960
P10141	Racked Micro Tips 2-200 µl (Bevelled Graduated) Natural, No. Per Case -10 Racks x 96 Pcs	960
P10107	Racked Micro Tips 2-200 µl (Xtreme Point Tip Orifice) Yellow, No. Per Case -10 Racks x 96 Pcs	960
P10108	Racked Micro Tips 2-200 µl (Xtreme Point Tip Orifice) Natural, No. Per Case -10 Racks x 96 Pcs	960
P10142	Racked Micro Tips 2-300 µl (Graduated) Natural, No. Per Case -10 Racks x 96 Pcs	960
P10109	Racked Micro Tips 100-1000 µl Blue, No. Per Case -10 Racks x 96 Pcs	960

P10110	Racked Micro Tips 100-1000 µl Natural, No. Per Case - 10 Racks x 96 Pcs	960
P10155	Racked Micro Tips 100-1000 µl Natural Natural, No. Per Case - 10 Racks x 96 Pcs	960
P10157	Racked Micro Tips 100-1250 µl (Graduated) Natural, No. Per Case - 10 Racks x 96 Pcs	960
P10159	Racked Macro Tips 100-5000 µl (Graduated) Natural, No. Per Case - 6 Racks x 24 Pcs	144
P10160	Racked Macro Tips 500-10000 µl (Graduated), Natural, No. Per Case - 6 Racks x 24 Pcs	144
P10117	Racked Micro Tips 0.2-10 µl Sterile (Graduated) Sterile Natural, No. Per Case - 10 Racks x 96 Pcs	960
P10147	Racked Micro Tips 0.2-10 µl Extended (Graduated) Sterile Natural, No. Per Case - 10 Racks x 96 Pcs	960
P10148	Racked Micro Tips 2-200 µl (Bevelled Graduated) Sterile Natural, No. Per Case - 10 Racks x 96 Pcs	960
P10149	Racked Micro Tips 2-200 µl (Bevelled Graduated) Sterile Yellow, No. Per Case - 10 Racks x 96 Pcs	960
P10111	Racked Micro Tips 2-200 µl (Xtreme Point Tip Orifice) Yellow Sterile Yellow, No. Per Case - 10 Racks x 96 Pcs	960
P10112	Racked Micro Tips 2-200 µl (Xtreme Point Tip Orifice) Natural Sterile Natural, No. Per Case - 10 Racks x 96 Pcs	960
P10150	Racked Micro Tips 2-300 µl (Graduated) Natural Sterile Natural, No. Per Case - 10 Racks x 96 Pcs	960
P10113	Racked Micro Tips 100-1000 µl Sterile Blue Sterile Blue, No. Per Case - 10 Racks x 96 Pcs	960
P10114	Racked Micro Tips 100-1000 µl Sterile Natural Sterile Natural, No. Per Case - 10 Racks x 96 Pcs	960
P10156	Racked Micro Tips 100-1000 µl XL (Graduated), No. Per Case - 10 Racks x 96 Pcs	960
P10158	Racked Micro Tips 100-1250 µl (Graduated), No. Per Case - 10 Racks x 96 Pcs	960
P10171	Racked Macro Tips 100-5000 µl, No. Per Case - 6 Racks x 96 Pcs	576
P10178	Racked Macro Tips 500-10000 µl (Graduated), Natural, No. Per Case - 6 Racks x 96 Pcs	576
P10118	Refill Micro Tips 0.2-10 µl (Graduated) Natural, No. Per Case - 12 Refill Racks x 96 Pcs	1152
P10152	Refill Micro Tips 2-200 µl (Bevelled Graduated) Natural, No. Per Case - 10 Refill Racks x 96 Pcs	960
P10153	Refill Micro Tips 2-200 µl (Bevelled Graduated) Yellow, No. Per Case - 10 Refill Racks x 96 Pcs	960
P10119	Refill Micro Tips 2-200 µl (Xtreme Point Tip Orifice) Yellow, No. Per Case - 10 Refill Racks x 96 Pcs	960
P10120	Refill Micro Tips 2-200 µl (Xtreme Point Tip Orifice) Natural, No. Per Case - 10 Refill Racks x 96 Pcs	960
P10154	Refill Micro Tips 2-300 µl (Graduated) Natural, No. Per Case - 10 Refill Racks x 96 Pcs	960
P10176	Refill Micro Tips 100-1000 XL µl (Graduated) Natural, No. Per Case - 6 Refill Racks x 96 Pcs	576
P10177	Refill Micro Tips 100-1250 µl (Graduated) Natural, No. Per Case - 6 Refill Racks x 96 Pcs	960

What Abdos Life Sciences stands for?



Precision



Reliability



Responsiveness



Innovation



Sustainability



Relationships

Get in touch with us:

 **011-41081897**

 **labtech@abdosindia.com**

 **www.abdoslifesciences.com**