

# Technical Data Sheet

Doc No–Ab/TDS/UTL-226

Rev No: 03 Date 11/2019

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Abdos Brand – Derma Gaurd White Latex Gloves – 9.5 inches length**

**Valid for Catalog No. U20344, U20345, U20346, U20449, U20450**

### **1. Description/Specification**

**1.1 Description:** Abdos Brand – Derma Guard White Latex Gloves – 9.5 inches length

1.1.1 These gloves are ideal product to protect you exposure to bacteria, viruses and chemicals. It also enhances your ability to offer maximum security against cross contamination.

1.1.2 Ambidextrous, textured fingertips and beaded cuff enhances wet & dry grip.

1.1.3 Softness provides superior comfort & natural fit.

1.1.4 These gloves are powder free, minimizing the potential for powder related complications, irritation and dermatitis.

1.1.5 Lower level of extractable proteins and chemical residue, hence lower incidences of allergic reactions.

1.1.6 Conforms to ASTM D6319 and EN 455 Standards.

**1.2 Design:** Hand gloves design.

### **1.3 Performance Specifications:**

<b>Code</b>	<b>Size</b>	<b>No. Per Case</b>
<b>U20344</b>	Small	100
<b>U20345</b>	Medium	100
<b>U20346</b>	Large	100
<b>U20449</b>	X Small	100
<b>U20450</b>	X Large	100

**1.4 Material:** Latex gloves



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 & 13485 Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903  
**GST :** 05AAGCA8605A1Z7

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com

# Technical Data Sheet

**1.5 Color:** White

**1.6 Quality control:**

Testing of various characteristics in accordance with the valid Specification


**1.7 Other information**

1.7.1 These gloves are used for working on lab, machine and any other work.

1.7.2 These gloves are ideal product to protect you any of microbial contamination.

1.7.3 These Flexible gloves are fit in hand & easy to use on working.

1.7.4 Powder free Chlorinated gloves



**ABDOS**

**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 & 13485 Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903  
**GST :** 05AAGCA8605A1Z7

**P :** +91 (0) 1332 235277  
**E :** labtech@abdosindia.com  
**W:** www.abdoslifesciences.com