

Technical Data Sheet

Doc No – Ab/TDS/LH-39

Rev No: 00 Date 02/2022-23

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Elegant Single Channel Adjustable Volume Pipettes, PP

Valid for Catalog No. E11150, E11151, E11152, E11153, E11154, E11155, E11156, E11157, E11158

1. Description/Specification

1.1 Description: Elegant Single Channel Adjustable Volume Pipettes, PP

1.1.1 Abdos Elegant pipettes are very accurate and gives desired results.

1.1.2 Abdos Elegant pipettes are lightweight, ergonomic, with low force design.

1.1.3 Volume settings can be clearly read through digital display with 4 digit counter.

1.1.4 Abdos Elegant pipettes are easy to calibrate & maintain with the tool supplied.

1.1.5 Abdos Elegant pipettes are designed to avoid repetitive strain injuries.

1.1.6 Abdos Elegant pipettes are supplied with operator manual that provides complete information on calibration & maintenance of the pipettes.

1.1.7 Abdos Elegant pipettes are suitable for transfer liquid volumes.

1.1.8 Matching universal style pipette tips, calibrated in accordance with ISO 8655. The quality control involves gravimetric testing of each pipette with distilled water at 22°C. Abdos introduces fully autoclavable, ergonomically designed Elegant micropipettes which offers high level of accuracy and precision. All Elegant Pipettes have Chemical & UV Resistant Body.

1.2 Design:

Light weight & simple to use operation.

1.3 Performance Specifications:

Push Button, Cap, Push Button Rod, Tip ejector, Grip cover, Handle, Display cover, Tip ejector collar, Tip cone, Volume setting, Filter element, Sealing & ejecting tips

1.4 Material: PP (polypropylene)

1.5 Color: White



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

1.6 Basic Features

1.6.1 Made up of highly Virgin thermoplastic polypropylene

1.6.2 CE certified

1.6.3 **Chemical resistance:** see website (www.abdoslifesciences.com)

1.7. Pipettes

Code	Capacity	No. Per Case
E11150	0.1-2.5 µl	1
E11151	0.5-10 µl	1
E11152	2-20 µl	1
E11153	10-100 µl	1
E11154	20-200 µl	1
E11155	30-300 µl	1
E11156	100-1000 µl	1
E11157	1000-5000 µl	1
E11158	2000-10000 µl	1

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakshewari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email: labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet



Abdos Labtech Pvt. Ltd.
Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
Email : labtech@abdosindia.com, www.abdoslifesciences.com
Ph: +91] 11 41081897



Elegant

Fully Autoclavable Micropipettes

Abdos introduces fully autoclavable, ergonomically designed Elegant micropipettes which offers high level of accuracy and precision. All Elegant Pipettes have Chemical & UV Resistant Body.

“NEW ADDITION”

LIQUID HANDLING

One of The
World's
Lightest
Pipette



Volume Lock



Lightweight, Ergonomic,
Low force design.



4 Digit Counter



Filter Element



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakshwari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7