

Technical Data Sheet

Doc No– Ab/TDS/LH-72

Rev No: 02 Date 01/2021

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Fast Release Pipette Pump

Valid for Catalog No. E11219, E11220, E11221

1. Description/Specification

1.1 Description: Fast Release Pipette Pump

1.1.1 Fast Release Pipette Pump ensures quick & accurate operation.

1.1.2 These light weight pipette pump have a fast release trigger that dispenses the entire pipette content.

1.1.3 The thumb wheel easily rotates for precision aspirating or dispensing, down to a volume of one drop at a time.

1.1.4 The thumb wheel can be used to blow out the contents of the pipette.

1.1.5 Resistant to acids, alkalis & easily disassembly for cleaning.

1.2 Design:

Fast Release Pipette Pump

1.3 Performance Specifications:

Code	Capacity (ml)	No. per Case
E11219	2	4
E11220	10	4
E11221	25	4



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet

1.4 Quality control:


Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.5. Basic Features:

1.5.1 Autoclavability: NO, Non autoclavable

1.5.2 Chemical resistance: see website (www.abdoslifesciences.com)

1.5.3 CE certified



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903