

# Technical Data Sheet

Doc No–Ab/TDS/UTL-211

Rev No: 05 Date 02/2022-23

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Glove Box Holder, Acrylic**

**Valid for Catalog No. P70301, P70302, P70303**

### **1. Description/Specification**

**1.1 Description:** Glove Box Holder, Acrylic

1.1.1 They designate a place for gloves.

1.1.2 Clear and pristine holder to give clean appearance.

1.1.3 Top loading of the glove box.

1.1.4 Easy to dispense the gloves through the holder.

1.1.5 Each holder has holes that help it to be installed in a convenient location to the proximity of your work and task.

**1.2 Design:** Glove Box Holder

**1.3 Performance Specifications:**

<b>Code</b>	<b>No. of Places</b>	<b>L x W x H (mm)</b>	<b>No. Per Case</b>
<b>P70301</b>	1	260 x 100 x 138	1
<b>P70302</b>	2	260 x 100 x 270	1
<b>P70303</b>	3	260 x 100 x 415	1

**1.4 Material:** Acrylic

**1.5 Color:** Transparent Natural

**1.6 Quality control:**

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakhewari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

**Ph. No: +91- (0) 1332235277**  
**Email:labtech@abdosindia.com**  
**Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)**  
**GST: 05AAGCA8605A1Z7**

# Technical Data Sheet

## 1.7 Other information


- 1.7.1 Excellent clarity and visibility.
- 1.7.2 Strong and lightweight material.
- 1.7.3 Good impact strength.
- 1.7.4 **Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))
- 1.7.5 Manufactured in Aseptic Environment

## Instruction for Use

Product can be used for food contact also.  
Storage temp. Max 30°C  
Before using the product check crackness etc.  
Check all the information printed on pouches or boxes.  
Always wear gloves while handling specimen.

## WARNING

**DO NOT** store strong oxidizing agents.  
**DO NOT** place any plastic labware in a flame.  
**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure  
Consult local fire codes prior to storage of flammable liquids.



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakshewari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

**Ph. No: +91- (0) 1332235277**  
**Email:labtech@abdosindia.com**  
**Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)**  
**GST: 05AAGCA8605A1Z7**