

Technical Data Sheet

Doc No–Ab/TDS/GLP-186

Rev No: 03 Date 01/2021

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Microscope Slide File, HIPS

Valid for Catalog No. P90113

1. Description/Specification

1.1 Description: Microscope Slide File, HIPS

1.1.1 These Slide File will hold up to 20 standard microscope slides 75 x 25 mm.

1.1.2 These slide file is made of two parts the base holding the slide horizontally in pre-numbered spaces for easy identification and the transparent door which can both cover the slide and allow easy reading of ID labels with or without optical bar code reader.

1.1.3 Removal of the slides is made easy by pressing on one end that will automatically lift the other end.

1.2 Design:

Slide File design

1.3 Performance Specifications:

Code	Places	L x W x H (mm)	No. Per Case
P90113	20	192 x 295 x 11	10

1.4 Material:

1.4.1 HIPS (high Impact Polystyrene) confirm US FDA 21CFR

1.4.2 Complies with requirement accordance with article of regulation (EC) No. 1907/2006 and its amendments.

1.4.3 Complies with the requirement of European Pharmacopoeia, 7th edition (2011) & its supplements 7.5 (07/2012) monographs 3.2.2 Plastic Containers and closures for Pharmaceutical use.

1.4.4 The product has passed the USP testing including class VI tests & has been assigned the FDA Drug master file DMF 12123



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A127

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet

1.4.5 Materials has successful passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.

1.4.6 Free from Heavy metals, RoHS Free according to meets relevant requirements of EU regulations 1935/2004/CE, 2002/95/EC as amended, UE 10/2011 as amended. EU CONEG & it's following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water

1.4.7 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our Knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product

1.5 Color: White

1.6 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release amendments.

1.7 Other information

1.7.1 Impact resistant and tough.

1.7.2 Temperature range from -80°C to 80°C

1.7.3 Superior for its hardness, gloss, toughness, and electrical insulation properties.

1.8 Basic Features:

1.8.1 Autoclavability: NO

1.8.2 Chemical resistance: see website (www.abdoslifesciences.com)

1.8.3 Manufactured in Aseptic Environment




Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001, 13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A127

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India

An ISO 9001 ,13485, 14001 & CE Certified

A Division of Abdos

Regd. Office:

13A, Government Place (East)
Kolkata – 700069, India

GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277

E : labtech@abdosindia.com

W: www.abdoslifesciences.com

CIN : U85100WB2008PTC123903