

Technical Data Sheet

Doc No – Ab/TDS/LH-39

Rev No: 03 Date 10/2021-22

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Premium Plus 12 Channel Pipettes, PP

Valid for Catalog No. E11139, E11140, E11141

1. Description/Specification

1.1 Description: Premium Plus Multi-Channel Pipettes, PP

1.1.1 Abdos Premium Plus Multi-Channel Pipettes are fully autoclavable.

1.1.2 Abdos Premium Plus Multi-Channel Pipettes are lightweight, ergonomic, with low force design.

1.1.3 Volume settings can be clearly read through digital display.

1.1.4 Abdos Premium Plus Multi-Channel Pipettes are easy to calibrate & maintain with the tool supplied.

1.1.5 Abdos Premium Plus Multi-Channel Pipettes are designed to avoid repetitive strain injuries.

1.1.6 Abdos Premium Plus Multi-Channel Pipettes are supplied with operator manual that provides complete information on calibration & maintenance of the pipettes.

1.1.7 Abdos Premium Plus Multi-Channel Pipettes are suitable for transfer liquid volumes.

1.1.8 Matching universal style pipette tips, calibrated in accordance with ISO 8655/DIN 12650. The quality control involves gravimetric testing of each pipette with distilled water at 22°C.

1.2 Design:

Light weight & simple to use operation.

1.3 Performance Specifications:

Push Button, Cap, Push Button Rod, Tip ejector, Grip cover, Handle, Display cover, Tip ejector collar, Tip cone, Volume setting, Sealing & ejecting tips

1.4 Material: PP (polypropylene)



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet

1.5 Color: White

1.6 Basic Features

1.6.1 Made up of highly Virgin thermoplastic polypropylene

1.6.2 Available for standard 96 well plates.

1.6.3 Dispensing head of these pipettes rotate for optimum pipetting convenience.

1.6.4 The special compound material tip cone design allows visual seal verification and has excellent chemical resistance.

1.6.5 Individual piston and tip cone assemblies allow easy repair and maintenance.

1.6.6 12- Channel adjustable premium plus pipettes

1.6.7 Autoclavability: Yes


1.6.8 CE Certified

1.6.9 Chemical resistance: see website (www.abdoslifesciences.com)

1.7. Pipettes

Code	Capacity	No. Per Case
E11139	0.5-10 µl	1
E11140	5-50 µl	1
E11141	50-300 µl	1

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS 

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903