

Technical Data Sheet

Doc No–Ab/TDS/UTL-203

Rev No: 04 Date 10/2021-22

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Resealable Bags with Zip Lock, PE

Valid for Catalog No. U40110, U40111

1. Description/Specification

1.1 Description: Resealable Bags with Zip Lock, PE

1.1.1 These bags can be used to store and ship samples, solids and food.

1.1.2 Suitable for use with food.

1.1.3 Though it can be used for all applications it has especially used in food, water, forensic science and quality control laboratory.

1.1.4 Using the rim and channel closure the bag can be closed by figure pressure and reopened repeatedly.

1.1.5 Content stored in the bag would be dust free and free of contaminates.

1.2 Design:

Rectangular shaped Resealable Bags.

1.3 Performance Specifications:

Code	Size WxL (inch)	No. Per Case
U40110	9.8 x 10.2	100
U40111	13 x 13.35	100

1.4 Material:

1.4.1 PE (Polyethylene), Hydrophobic, confirm US FDA 21CFR

1.4.2 Complies with requirement accordance with article of regulation (EC) No. 1907/2006 and its amendments.



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet

1.4.3 Complies with the requirement of European Pharmacopoeia, 7th edition (2011) & its supplements 7.5 (07/2012) monographs 3.2.2 Plastic Containers and closures for Pharmaceutical use.

1.4.4 The product has passed the USP testing including class VI tests & has been assigned the FDA Drug master file DMF 12123

1.4.5 Materials has successful passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.

1.4.6 Free from Heavy metals, RoHS Free according to meets relevant requirements of EU regulations 1935/2004/CE, 2002/95/EC as amended, UE 10/2011 as amended. EU CONEG & it's following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water

1.4.7 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our Knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color: Transparent Natural

1.6 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7 Other information


1.7.1 Better chemical stability.

1.7.3 Highly transparent, elastic thermoplastic.

1.7.4 Autoclave - NO

1.8 Chemical resistance: See website (www.abdoslifesciences.com)

1.9 Manufactured in Aseptic Environment



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903