

Technical Data Sheet

Doc No–Ab/TDS/UTL-227

Rev No: 03 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Safeskin Purple/Blue Nitrile Gloves 9.5 & 12.0 inches length

Valid for Catalog No. U20301, U20301B, U20302, U20302B, U20303, U20303B, U20304, U20305, U20306

1. Description/Specification

1.1 Description: Safeskin Purple/Blue Nitrile Gloves

1.1.1 These gloves are ideal product to protect you exposure to bacteria, viruses and chemicals. It also enhances your ability to offer maximum security against cross contamination.

1.1.2 Ambidextrous, textured fingertips and beaded cuff enhances wet & dry grip.

1.1.3 Softness provides superior comfort & natural fit.

1.1.4 These gloves are powder free, minimizing the potential for powder related complications.

1.1.5 Safeskin Blue Nitrile Gloves Powder free- sample meets the overall migration requirement under Commission Regulation (EU) No 10/2011- Plastic materials and articles shall not transfer their constituents to foodstuff in quantities exceeding 10mg of total constituents released per dm² of food contact surface (mg/dm²) (overall migration limit)

1.2 Design: Hand gloves design

1.3 Performance Specifications:

1.3.1 Safe skin Purple Nitrile Gloves – 9.5 inches length

Code	Size	No. Per Case
U20301	Small	100
U20302	Medium	100
U20303	Large	100



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.3.1 Safe skin Purple Nitrile Gloves – 12.0 inches length

Code	Size	No. Per Case
U20304	Small	50
U20305	Medium	50
U20306	Large	50

1.3.3 Safe skin Blue Nitrile Gloves – 9.5 inches length

Code	Size	No. Per Case
U20301B	Small	100
U20302B	Medium	100
U20303B	Large	100

1.4 Material: Nitrile gloves

1.5 Color: Purple & Blue

1.6 Quality control:

Testing of various characteristics in accordance with the valid Specification


1.7 Other information


1.7.1 These gloves are used for working on lab, machine and any other work.

1.7.2 These gloves are ideal product to protect you any of microbial contamination.

1.7.3 These Flexible gloves are fit in hand & easy to use on working.

1.7.4 Powder free Chlorinated gloves


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com