

# Technical Data Sheet

Doc No–Ab/TDS/UTL-207

Rev No: 03 Date 01/2021

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Safety Face Shield, PC**

**Valid for Catalog No. U40102**

### **1. Description/Specification**

**1.1 Description:** Safety Face Shield, PC

1.1.1 Clear face shield of polycarbonate provides high impact and splash protection.

1.1.2 Impact resistant and lightweight frames for comfortable Face safety.

1.1.3 Ratched adjustment on the headgear provides comfortable fit.

**1.2 Design:** Safety Face Shield

**1.3 Performance Specifications:**

Code	L x B x T (inch)	No. Per Case
U40102	15.5 x 8 x 0.04	1

**1.4 Material:** PC (Polycarbonate)

**1.5 Color:** Top: Black      **Frame:** Transparent Natural

**1.6 Quality Control:**

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

**1.7 Other Information**

1.7.1 Excellent clarity and visibility.

1.7.2 Moderate chemical stability.

1.7.3 Transmits light better than alternative material.



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttarakhand, India  
**An ISO 9001 ,13485, 14001 & CE Certified**  
A Division of Abdos


**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**GST : 05AAGCA8605A1Z7**

**P : +91 (0) 1332 235277**  
**E : labtech@abdosindia.com**  
**W: www.abdoslifesciences.com**  
**CIN : U85100WB2008PTC123903**

# Technical Data Sheet

**1.8 Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.9 Manufactured in Aseptic Environment**



**ABDOS**

**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India

**An ISO 9001 ,13485, 14001 & CE Certified**

A Division of Abdos

**Regd. Office:**

13A, Government Place (East)  
Kolkata – 700069, India

**GST : 05AAGCA8605A1Z7**

**P : +91 (0) 1332 235277**

**E : labtech@abdosindia.com**

**W: www.abdoslifesciences.com**

**CIN : U85100WB2008PTC123903**