

Technical Data Sheet

Doc No–Ab/TDS/UTL-221

Rev No: 02 Date 11/2019

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Safety Goggles Box, Acrylic

Valid for Catalog No. P70106, P70107

1. Description/Specification

1.1 Description: Safety Goggles Box, Acrylic

1.1.1 This safety box can be set on the bench top or mounted on the wall for storage of goggles, glasses and other small lab product.

1.1.2 Safety Lid Keep the box free of dust and contamination.

1.1.3 Clear window helps to view the inventory from the outside.

1.2 Design: Safety Goggles Box

1.3 Performance Specifications:

Code	L x W x H (cm)	Size	No. Per Case
P70106	23 x 16 x 18	Small	1
P70107	23 x 16 x 25	Large	1

1.4 Material: Acrylic

1.5 Color: Transparent Natural

1.6 Quality control:

Testing of various characteristics in accordance with the valid Specification



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos


Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com

Technical Data Sheet

1.7 Other information

- 1.7.1 Excellent clarity and visibility.
- 1.7.2 Strong and lightweight material.
- 1.7.3 Good impact strength.



ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 & 13485 Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com