

**ABDOS**



# Serological Pipettes

---

Creating For Life Sciences

[www.abdoslifesciences.com](http://www.abdoslifesciences.com)



- Sterile Serological has been validated according to EN ISO 11137-2; 2012 for SAL  $10^{-6}$ .
- Advance PU filter offers complete protection against contamination.
- Bidirectional graduations provide choice of ascending and descending scales.
- Color-coding on paper/serological pipettes allows easy size identification.
- Negative graduations provide additional pipetting volumes.
- Accuracy within  $\pm 2\%$  at full volume.
- Individual lot number for traceability.



Cotton plug can be a major source of contamination to your experiment. Abdos advance PU filter offers complete protection against contamination which gives hassle & anxiety free results.

Code	Length (mm)	Capacity (ml)	No.Per Case
PI0601	262	1	500
PI0602	280	2	400
PI0603	333	5	200
PI0604	354	10	200
PI0605	312	25	200
PI0606	360	50	100

Abdos Serological Pipettes are tested for endotoxins (<0.005 EU/mL), Cytotoxicity (ISO 10993-5) and are guaranteed non-haemolytic.

For additional product or technical information,  
please visit [www.abdoslifesciences.com](http://www.abdoslifesciences.com).  
For distributor information contact [labtech@abdosindia.com](mailto:labtech@abdosindia.com).