

# Technical Data Sheet

Doc No–Ab/TDS/GLP-186

Rev No: 05 Date 02/2022-23

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Slide Box, ABS**

**Valid for Catalog No. P90103, P90104, P90105**

### **1. Description/Specification**

#### **1.1 Description:** Slide Box, ABS

1.1.1 These slide tray is suitable for holding standard 25 x 75 mm microscope slides.

1.1.2 Slide tray is mounted with cock sheet on the bottom and index list inside the lid for microscope slide identification.

1.1.3 Uniform grooved slide slots separates slides and prevents it to come in contact with each other.

1.1.4 Hinged lid gives secure closure to the box.

#### **1.2 Design:**

Rectangular shape Box

#### **1.3 Performance Specifications:**

| <b>Code</b>   | <b>Places</b> | <b>L x W x H (mm)</b> | <b>No. Per Case</b> |
|---------------|---------------|-----------------------|---------------------|
| <b>P90103</b> | 25            | 95 x 82 x 32          | 12                  |
| <b>P90104</b> | 50            | 199 x 82 x 32         | 6                   |
| <b>P90105</b> | 100           | 192 x 162 x 32        | 6                   |

#### **1.4 Material:**

1.4.1 ABS (Acrylonitrile Butadiene Styrene) confirm US FDA 21CFR

1.4.2 Complies with the requirement of European Pharmacopoeia, Plastic Containers and closures for Pharmaceutical use.

1.4.3 The product has passed the USP testing including class VI tests ,Materials has successful passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakheswari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

**Ph. No: +91- (0) 1332235277**  
**Email:labtech@abdosindia.com**  
**Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)**  
**GST: 05AAGCA8605A1Z7**

# Technical Data Sheet

1.4.4 Free from Heavy metals, RoHS Free according to meets relevant requirements of EU regulations 1935/2004/CE, 2002/95/EC as amended, UE 10/2011 as amended. EU CONEG & its following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water

1.4.5 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our Knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

**1.5 Color:** Grey

## **1.6 Quality control:**

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release amendments.

## **1.7 Other information**

1.7.1 Impact resistant and tough.

1.7.2 Temperature range from -20°C to 80°C

1.7.3 Superior for its hardness, gloss, toughness, and electrical insulation properties.

## **1.8 Basic Features:**

**1.8.1 Autoclavability:** NO

**1.8.2 Chemical resistance:** see website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.8.3** Manufactured in Aseptic Environment

## **Instruction for Use**

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

## **WARNING**

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.

**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



Abdos Labtech Private Limited  
(A Division of Abdos Group)  
(An ISO 9001, 13485, 14001, 18385, GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakheswari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277  
Email: labtech@abdosindia.com  
Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)  
GST: 05AAGCA8605A1Z7

# Technical Data Sheet



**Abdos Labtech Pvt. Ltd.**

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India  
Email: labtech@abdosindia.com, www.abdoslifesciences.com  
Ph: +91 11 41081897



Quality Assurance Incharge  
Abdos Labtech Pvt. Ltd.



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakshwari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903

**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277  
Email:labtech@abdosindia.com  
Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)  
GST: 05AAGCA8605A1Z7