

# Technical Data Sheet

Doc No–Ab/TDS/UTL-227

Rev No: 05 Date 02/2022-23

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Specialty Indicator Tape for Steam Autoclave**

**Valid for Catalog No. U20802**

### **1. Description/Specification**

**1.1 Description:** Specialty Indicator Tape for Steam Autoclave

1.1.1 These autoclavable tapes adhere to almost any dry surface-glass, plastic, metal, paper and rubber.

1.1.2 Versatile tape indicates completion of normal steam sterilization process.

1.1.3 All specialty steam indicator tape accepts and written notations may be used to seal and identify almost all laboratory material before and after processing.

1.1.4 Change in the colour of the written matter from white to dark completion of autoclave cycle.

**1.2 Design:** Specialty Indicator Tape for Steam Autoclave

**1.3 Performance Specifications:**

<b>Code</b>	<b>Size (Inch)</b>	<b>No. Per Case</b>
<b>U20802</b>	1 x 500	1 Roll

**1.4 Color:** White

**1.5 Quality control:**

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

**1.6 Other information:**

1.6.1 **Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com)) 1.6.2 Manufactured in Aseptic Environment

### Instruction for Use

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakheswari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903

**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

**Ph. No: +91- (0) 1332235277**  
**Email:labtech@abdosindia.com**  
**Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)**  
**GST: 05AAGCA8605A1Z7**

# Technical Data Sheet

## WARNING

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.

**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

