

Technical Data Sheet

Doc No – Ab/TDS/Eqpt -232

Rev No : 06 Date 10/2021-22

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Swirltop – Ultra Flat & Eco Magnetic Stirrer

Valid for Catalog No. E11307, E11308

1. Description/Specification

1.1 Description: Swirltop – Ultra Flat & Eco Magnetic Stirrer

E11307 Ultra Flat Magnetic Stirrer

E11308 Eco Magnetic Stirrer

1.1.1 SWIRLTOP magnetic stirrer is used in biotech labs, provides better mixing & is ideal for application of physical & chemical analysis.

1.1.2 Its body is made of ABS material.

1.1.3 It comes with automatic reverse rotation for better mixing result.

1.1.4 Top plate is made from chemically resistant materials (PET).

1.2 Design: Swirltop – Ultra Flat & Eco Magnetic Stirrer

1.3 Performance Specifications:

Code	E11307	E11308
Voltage	220-240V	220-240V
Frequency	50-60Hz	50-60Hz
Motor Type	Motor less	Brushless DC Motor
Stirring position	01	01
Speed range	15–1500 rpm	300–2000 rpm
Max. stirring quantity (H ₂ O)	0.5 Ltr	1.5 Ltr
Maximum magnetic bar (length)	40mm	40mm
Dimensions of work plate	90x90mm	120x120mm
Dimensions	175x120x15 mm	150x145x50 mm
Weight	0.3 kg	0.3 kg
Speed setting	Button	Scale



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet

Permissible ambient Temperature	5-40°C	5-40°C
Permissible relative humidity	80%	80%
Protection Class	P65	P21

1.4 Material: Body – ABS Material

1.5 Color: Grey

1.6. Basic features:

1.6.1 **Material:** Made up of Highly pure Virgin ABS (Acrylonitrile Butadiene Styrene)

1.6.2 CE certified

1.6.3 **Chemical resistance:** see website (www.abdoslifesciences.com)

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.




ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903