

Technical Data Sheet

Doc No – Ab/TDS/Eqpt -230

Rev No : 07 Date 02/2022-23

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Swirltop – Ultra Flat & Eco Magnetic Stirrer

Valid for Catalog No. E11307, E11308, E11309

1. Description/Specification

1.1 Description: Swirltop – Ultra Flat & Eco Magnetic Stirrer

E11307 Ultra Flat Magnetic Stirrer

E11308 Eco Magnetic Stirrer Round

E11309 Eco Magnetic Stirrer Aquare

1.1.1 SWIRLTOP magnetic stirrer is used in biotech labs, provides better mixing & is ideal for application of physical & chemical analysis.

1.1.2 Its body is made of ABS material.

1.1.3 It comes with automatic reverse rotation for better mixing result.

1.1.4 Top plate is made from chemically resistant materials (PET).

1.2 Design: Swirltop – Ultra Flat & Eco Magnetic Stirrer

1.3 Performance Specifications:

Code	E11307	E11308	E11309
Voltage	220-240V	220-240V	220-240V
Frequency	50-60Hz	50-60Hz	50-60Hz
Motor Type	Motor less	Brushless DC Motor	Brushless DC Motor
Stirring position	01	01	01
Speed range	15–1500 rpm	300–2000 rpm	300–2000 rpm
Max. stirring quantity (H2O)	0.8 Ltr	1.5 Ltr	1.5 Ltr
Maximum magnetic bar (length)	40mm	40mm	40mm
Dimensions of work plate	90x90mm	120x120mm	120x115mm
Dimensions	175x120x15 mm	150x145x50 mm	150x144x50 mm
Weight	0.3 kg	0.3 kg	0.3 kg
Speed setting	Button	Scale	Scale
Permissible ambient Temperature	5-40°C	5-40°C	5-40°C



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email: labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

Permissible relative humidity	80%	80%	80%
Protection Class	IP65	IP21	IP21

1.4 Material: Body – ABS Material

1.5 Color: Grey

1.6. Basic features:

1.6.1 Material: Made up of Highly pure Virgin ABS (Acrylonitrile Butadiene Styrene)

1.6.2 CE certified

1.6.3 Chemical resistance: see website (www.abdoslifesciences.com)

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing

Instruction for Use

Storage temp. Max 30°C

Before using the product check crackness, etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



Abdos Labtech Pvt. Ltd.

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
Email : labtech@abdosindia.com, www.abdoslifesciences.com
Ph: +(91) 11 41081897



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakshwari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslifesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email: labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7