

Technical Data Sheet

Doc No Ab/TDS/Eqpt-251

Rev No: 04 Date 02/2022-23

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

WAVEMAX – Digital Incubator

Valid for Catalog No. E11355, E11356

1. Description Specifications:

1.1 Description: WAVEMAX – Digital Incubator

1.1.1 Wavemax-Digital Incubators provide the facility to monitor the chamber temperatures in real time.

1.1.2 Has large internal chambers with two adjustable shelves for increased space. Internal Power Outlet for Installation of other

1.1.3 Equipment for sample agitation within the Digital Incubators.

1.1.5 Wavemax-Digital Incubator is available in a heating unit or heating and cooling for lower temperature applications.

1.2 Design: Digital Incubator

1.3 performance Specifications:

Specifications	E11355	E11356
Voltage	220-240V	220-240V
Frequency	50-60 Hz	50-60 Hz
Temp. range	Ambient +5 to 60°C	Ambient +5 to 60°C
Platform Dimensions	235x290mm	235x290mm
Operating Temp. range	+4 to 65°C	+4 to 65°C
Exterior	335x370x475mm	335x370x475mm
Interior	26x23.5x32.5 cm	26x23.5x32.5 cm
Capacity	20L	20L
Weight	6.5kg	6.5kg

1.4 Material: Body – ABS

1.5 Color: Grey

1.6 E11355 – Wavemax – Heating Digital Incubator

1.7 E11355 – Wavemax – Heating & Cooling Digital Incubator



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7

Technical Data Sheet

1.8 Spare Accessories

1.8.1 **E11357** – Extra shelf 10.5x8 inch plastic (not autoclavable)

1.8.2. CE certified

1.8.3 **Chemical resistance:** see website (www.abdoslifesciences.com)

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing.

Instruction for Use

Storage temp. Max 30°C

Before using the product check crackness, etc.

Check all the information printed on pouches or boxes.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



Technical Data Sheet



Abdos Labtech Private Limited
(A Division of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903

Corporate Office
319, 3rd Floor, DLF Tower A, DDA Dist. Centre
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277
Email:labtech@abdosindia.com
Web: www.abdoslivesciences.com
GST: 05AAGCA8605A1Z7