

Technical Data Sheet

Doc No Ab/TDS/Eqpt-244

Rev No : 04 Date 10/2021-22

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

WAVEX – CLASSIC TUBE ROTATOR

Valid for Catalog No. E11300

1. Description/Specification

1.1 Description:

1.1.1 Wavex tube rotators provide gentle but effective mixing of biological samples, held in 1.5 ml to 50 ml tubes.

1.1.2 Used in a variety of applications including immune precipitations, prevention of blood coagulation & latex diagnostics.

1.1.3 Wavex-Classic Tube Rotator has Disk Type Design.

1.2 Design:

Disk type design

1.3 performance Specifications:

1.3.1 Voltage (AC)	: 220~240 Volt
1.3.2 Frequency (Hz)	: 50
1.3.3 Power (W)	: 20
1.3.4 Tilt Angle	: 0- 90°
1.3.5 Motor type	: DC Motor
1.3.6 Speed range (RPM)	: 0 ~ 80
1.3.7 Operation time	: Continuous
1.3.8 Dimensions	: 300x220x310mm
1.3.9 Permissible ambient Temp	: 5- 40°
1.4.0 Permissible Relative humidity:	80 %
1.4.1 Speed Display :	Scale



Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttrakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903

Technical Data Sheet

1.4 Material: Body – ABS, Stainless Steel

1.5 Color: Dark Blue


1.6. Spare Accessories:

E11301 Plastic clips for 1.5ml x 12 MCT, Disk accessory – 12pcs/case
E11302 Plastic clips for 15ml x 12 MCT, Disk accessory – 12pcs/case
E11303 Plastic clips for 50ml x 12 MCT, Disk accessory – 12pcs/case

1.7. CE certified

1.8 **Chemical resistance:** see website (www.abdoslifesciences.com)

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.




ABDOS

Shiv Ganga Industrial Estate,
Plot no.1 Lakeshwari, Bhagwanpur,
Roorkee – 247667 Uttarakhand, India
An ISO 9001 ,13485, 14001 & CE Certified
A Division of Abdos

Regd. Office:
13A, Government Place (East)
Kolkata – 700069, India
GST : 05AAGCA8605A1Z7

P : +91 (0) 1332 235277
E : labtech@abdosindia.com
W: www.abdoslifesciences.com
CIN : U85100WB2008PTC123903